
1. Use Case

As a [Pathologist], I want to raise the re-stain request for the new slides preparation

2. Visual UI



3. Pre-condition

- ❏ Only the Result entry pending list for result entry should be visible.
- ❏ All the result entry list should be visible based on the logged unit.
- ❏ After acceptance all the tests should remain on this screen until the report will not be marked as final released.
- ❏ All the tests should be displayed based on the user sub department access role.
- ❏ Users can do multiple time data updates on this screen until the report release is not done.
- ❏ Users can only access those tests which are assigned to logged in units.

4. Acceptance Criteria

- ❏ Pathologist should have the facility to initiate the request for re-stain in some cases where pathologists want to reconfirm specimen reporting
- ❏ To initiate re-stain at the time of result reporting pathologist click on the restaining/reblocking option.

- Once user click on the restaining/reblocking button user will get 2 options i.e. Against Block & Against Slides.
- To initiate the restain user will select the Agaist Slide with Blocks :
  - **Default View**




- If user want to initiate the new slide request against the blocks then user will follow the below steps :
- At the block level once a user clicks on the [+] button the new line will be added under the same block for the newly requested slides with the following running number and default staining type.




- Once the new request is added before the save user has the option to remove.
- Users will also have the option to change the staining type if required.
- Once user clicks on the save application will send the new slides to the “Staining” to prepare the slides with the type = “ReStain”

- To initiate the restain user will select the Agaist Slide without Blocks :
  - **Default View**



  - If user want to initiate the new slide request without blocks then user will follow the below steps :
  - At the last slides level once a user clicks on the [+] button the new line will be added for the newly requested slides with the following running number and default staining type.



  - Once the new request is added before the save user has the option to remove.
  - Users will also have the option to change the staining type if required.
  - Once user clicks on the save application will send the new slides to the “Staining” to prepare the slides with the type = “ReStain”

5. Technical Requirement

- Application need to follow the running number for newly requested slides in the same running sequence with and without block
- Users will also have the facility to perform the re-staying for the new slides and existing slides together.

6. Outcome

- Once user clicks on the save application will send the new slides to the “Staining” to prepare the slides with the type = “ReStain”

As a [Pathologist], I want to raise the re-stain request for the new slides preparation

---

-: End of the Use Case :-