
1. Use Case

As a [Pathologist], I want to raise the re-stain request for the existing slide re-staining

2. Visual UI



3. Pre-condition

- ❑ Only the Result entry pending list for result entry should be visible.
- ❑ All the result entry list should be visible based on the logged unit.
- ❑ After acceptance all the tests should remain on this screen until the report will not be marked as final released.
- ❑ All the tests should be displayed based on the user sub department access role.
- ❑ Users can do multiple time data updates on this screen until the report release is not done.
- ❑ Users can only access those tests which are assigned to logged in units.

4. Acceptance Criteria

- ❑ Pathologist should have the facility to initiate the request for re-stain against the existing slides in some cases where pathologists want to reconfirm specimen reporting
- ❑ To initiate re-stain against the existing slides at the time of result reporting pathologist click on the restaining/reblocking option.

- Once user click on the restaining/reblocking button user will get 2 options i.e. Against Block & Against Slides.
- To initiate the restain against existing slided user will select the Agaist Slide with Blocks :
  - **Default View**



  - If user want to initiate the re-stain request for the existing slide new slide request against the slides then user will follow the below steps :
  - At the slides level once a user clicks on the [R] button the new line will be added with the same slide number and the user will get the option to select the required staining type.



  - Once the new request is added before the save user has the option to remove.
  - Users will also have the option to change the staining type if required.
  - Once user clicks on the save application will send the new slides to the “Staining” to prepare the slides with the type = “ReStain”
- Same process will be followed for the slides which are without any blocks and will follow the same process.

5. Technical Requirement

- ❏ Application need to follow the running number for newly requested slides in the same running sequence with and without block
- ❏ Users will also have the facility to perform the re-staying for the new slides and existing slides together.

6. Outcome

- ❏ Once user clicks on the save application will send the new slides to the “Staining” to prepare the slides with the type = “ReStain”

-: End of the Use Case :-